
Force-Displacement Glenoid Reamer Guide



Western  Engineering



Author:

MATTHEW STOKES

March 7, 2014

Contents

1	Overview	5
2	Equipment	5
2.1	Sensing and Tools	5
2.2	3D Printed Attachment Parts	5
2.3	DAQ & Other Intermediary Electronics	6
2.4	Software	6
2.5	Supplementary Software Files	6
2.6	Misc Cables, Nuts, Bolts, Wires	7
3	Hardware Setup	7
3.1	Optical System Setup	8
3.1.1	Setup Test	10
3.1.2	Offset Calibration	11
3.2	Load Cell Setup	15
3.2.1	Setup Test	16
4	Sensor Software	18
4.1	Return of the Motion Station	18
4.2	Loadcelly	21
4.3	Recording	22
5	Experimental Procedure	23
6	Post Processing	24
Appendix A Supplementary Software Files		27
A.1	smart_02	27
A.2	bsty0213.rig	28
A.3	FT13036Nano25ReCalibrated.cal	29
Appendix B Code		31
B.1	Packets	31
B.1.1	Packet.py	31
B.1.2	OpticalTrackerData.py	32
B.1.3	LoadCellData.py	34
B.2	Sensors	35
B.2.1	Sensor.py	35
B.2.2	LoadCell.py	36
B.2.3	OpticalTracker.py	38
B.3	Visitor Design Pattern	40
B.3.1	Reamer.py	40
B.3.2	DataExchange.py	42
B.3.3	Chart.py	43

	B.3.4	StiffnessElement.py	52
B.4	GUI	53
	B.4.1	gui.py	53
	B.4.2	main_gui.py	56

List of Figures

1	Reamer Outfitted with Load Cell and Optical Tracker	7
2	Optical Tracker System Overview	8
3	Optotrak Box	9
4	Certus Optical Camera	9
5	NDI First Principals	10
6	NDI First Principals	10
7	Experiment Setup	11
8	Marker Setup	12
9	Advanced Marker Settings	13
10	Rigid Body Setup Tab	13
11	Two Rigid Bodies	13
12	Wireless Mode	14
13	Offset	14
14	Load Cell Overview	15
15	ATI DAQ	15
16	Broken 24 Pin Cable	15
17	NI USB-6211	16
18	National Instruments Measurement & Automation Explorer	16
19	DEV port	17
20	ATIDAQFT.NET	17
21	DAQ Hardware Options	18
22	Return to the Motion Station	18
23	Create New Session Folder	19
24	Return to Motion Station Running	19
25	Add/Configure Tracker	19
26	Certus Setup	20
27	Calibration Matrix	21
28	Loadcelly	21
29	Loadcelly Calibration File	22
30	Recording	23
31	GUI	24
32	GUI Normalized	25
33	GUI Normalized and Biased	25
34	GUI Normalized, Biased, and Trimmed	26

List of Tables

1	Sensing and Tools	5
2	3D Printed Attachment Parts	5
3	DAQ & Other Intermediary Electronics	6
4	Software	6
5	Supplementary Software Files	6
6	Misc Cables, Nuts, Bolts, Wires	7
7	Reamer Assembly Parts List	8
8	Reamer Assembly Parts List	12
9	ot.csv	22
10	lc.csv	23
11	Excel	26

1 Overview

Reamer project was designed to calculate the applied force and resulting displacement of the glenoid during reaming. However due to modern sensor limitations (namely a small wireless load cell) the reamer will not be rotating during measurements.

This document includes information pertaining to setting up the system, recording measurements from the sensors, and interacting with this information. It does not include any background information about the project.

2 Equipment

Equipment has been separated into six distinctly different categories: sensing and tools, 3D printed attachment parts, DAQ & other intermediary electronics, software, supplementary files, and misc.

2.1 Sensing and Tools

Item	Product	Manufacture	Description
1	Load cell	ATI	Nano25 FT 13036
2	Optotrack Certus	NDI	Optical Camera
3	Optical tracker	NDI	
4	Reamer		Surgical reamer

Table 1: Sensing and Tools

2.2 3D Printed Attachment Parts

All parts were designed by Matthew Stokes and printed on his MakerBot Replicator using the settings 100% infill, 0.2mm layer height.

Item	Product	Description
5	Head Couple	Between head of reamer and load cell
6	Back Plate Couple	Between back plate of load cell and reamer shaft
7	Reamer Attachment Plate	Base connecting to the reamer
8	Optical Tracker Attachment Plate	connecting to the reamer attachment plate
9	Standoff	Hidden via thru hole securing reamer attachment plate

Table 2: 3D Printed Attachment Parts

2.3 DAQ & Other Intermediary Electronics

Necessary intermediary products connecting to sensors and involved in setup/calibration.

Item	Product	Manufacture	Description
10	Wireless Strober	NDI	Strobe the optical trackers
11	Stylist	NDI	tracker offset to tip of stylist
12	Optotrak Box	NDI	Sync strober and tracking camera
13	Load Cell DAQ	ATI	filter and amplify
14	Multifunction DAQ	NI	load cell DAQ to USB
15	Computer		Windows XP

Table 3: DAQ & Other Intermediary Electronics

2.4 Software

Item	Program	Author	Description
16	Motion Station	HULC Lab	LabView optical tracker measurements
17	Loadcelly	Stokes	LabView load cell measurements
18	pushpush.py	Stokes	Python interacting with data output
19	ATIDAQFT.NET	ATI	Native software
20	NDI First Principals	NDI	Native software
21	NI Measurement & Automation	NI	Native software

Table 4: Software

2.5 Supplementary Software Files

These files should not be modified.

Item	File Name	Extension	Description
22	nano25calibration	.cal	load cell calibration
23	bsty0213	.rig	rigid body for the stylist
24	smart_02	.rig	rigid body for base NI optical tracker

Table 5: Supplementary Software Files

2.6 Misc Cables, Nuts, Bolts, Wires

Item	Product	Description
25	Threaded Rob	Connecting head couple to load cell
26	4x Bolts	Connecting reamer attachment plate and optical tracker plate
27	4x Nuts	
28	2x Screws	Attach reamer attachment plate to standoff
29	Broken out 24 Pin Cable	Connect Load Cell DAQ to NI-USB-6211
30	Thread	Secure reamer attachment plate to back end of reamer

Table 6: Misc Cables, Nuts, Bolts, Wires

3 Hardware Setup

The project consists of two main systems which are the load cell and the optical tracker. Setup for each of these systems is independent. The reamer itself will be setup as shown in the image below.



Figure 1: Reamer Outfitted with Load Cell and Optical Tracker

Contained in this assembly are the following items:

Item	Product
1	Load Cell
3	Optical Tracker
4	Reamer
5	Head Couple
6	Back Plate Couple
7	Reamer Attachment Plate
8	Optical Tracker Plate
9	Standoff
25	Threaded Rod
26	4X Bolts
27	4X Nuts
28	2X Screws
30	Thread

Table 7: Reamer Assembly Parts List

3.1 Optical System Setup

An overview of the optical system setup is shown in the image below.




Figure 2: Optical Tracker System Overview

Plug in and power on the optical tracker camera. You should see a green light which indicates the power is on. The optical tracker camera communicates with the Optotrak Box which is connected via Ethernet to the main computer. Power on the Optotrak Box and plug the cable into the optical tracker.



Figure 3: Optotrak Box



Figure 4: Certus Optical Camera

Next we need to wire the tracker on our reamer. The tracker should connect to the strober box which can hold up to 5 trackers. Trackers may also be daisy chained with up to 4 trackers. To begin we will wire the tracker as part of our reamer assembly and also a stylist (which we be used for offset calibration). Note we will wire the stylist in port 1 and the tracker in port 2. This order matters later.

3.1.1 Setup Test

To ensure the system is setup properly open NDI First Principals and click *Create New*




Figure 5: NDI First Principals

If everything is connected we should will be taken to the next screen for setting up a new project seen below. Else go back and try again.




Figure 6: NDI First Principals

3.1.2 Offset Calibration

Now that we know our optical tracker is being read by the camera we will calculate the offset between the center of the tracker and the tip of the reamer. We do because we are interested in the displacement of the tip as opposed to the displacement of the tracker. This allows us to also pivot about the tip without measuring any displacement.

To determine the offset we will position the stylist to touch the tip of the reamer and record the x,y,z values of both the stylist and the other tracker. The difference between these values will be our offset. We will later offset the trackers calibration matrix by these values in return to motion station.

Clicking the *Next* button on the multiple position sensor registration screen should take us to the experiment setup screen shown below.




Figure 7: Experiment Setup

Here we should see **6 total markers**- 3 from our stylist and 3 from our tracker. If you do not see six total then something has not been setup correctly.

At this stage we are required to modify some of the settings in order to get better functionality out of our optical system. Modify the following parameters:

Parameter	New Value
Collection Frame Frequency	50
Frequency	3500
Advanced — Duty Cycle	80

Table 8: Reamer Assembly Parts List

These modifications are also shown in the images below.

The screenshot shows the 'Marker Setup' tab of a software interface. At the top, there are three tabs: 'Marker Setup', 'Rigid Body Setup', and 'ODAU Setup'. Below the tabs, a text box instructs the user to describe how individual markers are connected to the system, reminding them to account for empty marker channels on any connected Optotrak 3020 strober. Below this, there is a section titled 'Total Number of Markers Per Strober Port' containing three input fields: 'Port 1:' with the value '6', 'Port 2:' with the value '0', and 'Port 3:' with the value '0'. A 'Configure Tools/Strobers' button is located below these fields. Further down, the 'Collection Frame Frequency' is set to '50 Hz'. Below that, the 'Marker Power/Frequency' section features a slider for 'Power' ranging from 'min' to 'max', with a current value of '39 %'. The 'Frequency' is set to '3500 Hz'. An 'Advanced Marker Settings' button is located at the bottom right of this section.

Figure 8: Marker Setup

Once these settings have been changed we need to click on the Rigid Body Setup Tab.

Here we need to add two rigid body files which can be found in Appendix A: **bsty0213.rig** and **smart_02.rig**. Bsty2013 was created in house and smart_02 is the rigid body file for a tracker created by NDI. These files contain information used by the camera to interoperate 3 markers as a rigid body. You must ensure to load the bsty0213 file before the smart_02 file as this is order




Figure 9: Advanced Marker Settings




Figure 10: Rigid Body Setup Tab

corresponds to our placement earlier. The completed rigid body setup is shown below.




Figure 11: Two Rigid Bodies

We can now click *Next* as we are finished the experimental setup. Clicking *Next* may or may not

prompt you with the screen seen below. Click *no*. Just do it. You are lying to yourself if you believe this will actually work in wireless mode.




Figure 12: Wireless Mode

Green means the camera can see the markets and I assume you can figure out what red means. Interesting note: if you are in an odd situation where you think it should be reading values but it's not simply put you hand in front of the tracker to hide it from the camera then remove it. This does something along some lines of resetting some internal stuff.

Click the *View — Probe View*. This should present you with a screen similar to the one shown below (minus any actual values we're getting there).




Figure 13: Offset

This is where we determine the offset. We use the smart.02 sensor as the reference and the bsty0213 as the probe. Then once both rigid bodies are green (in sight) it should start showing values which are the difference between the two trackers. It's important to note that the stylist has already been configured to offset its tracker to the tip so measurements are from the tip of the stylist to the center of the optical tracker on the reamer.

Touch the tip of the stylist to the tip of the reamer and hit *F5* or the *Capture* button. The resulting values should be similar to the ones shown in the image above else something has really changed which is odd. Record these values somewhere as we will use these soon. It's at this time we are ready to move on to using Return of the Motion Station (yay). But before we do this we will take a second to setup our load cell.

3.2 Load Cell Setup

An overview of the optical system setup is shown in the image below.




Figure 14: Load Cell Overview

The load cell plugs into the ATI DAQ box. This box connects to the NI USB-6211 device via a special broken out cable.



Figure 15: ATI DAQ



Figure 16: Broken 24 Pin Cable

Finally, the NI USB-6211 box connects to the USB port on the main computer.



Figure 17: NI USB-6211

3.2.1 Setup Test

Open NI Measurement & Automation Explorer which should display the screen shown below.




Figure 18: National Instruments Measurement & Automation Explorer

Select *Devices and Interfaces* — *NI-DAQmx Devices* and make note of the **DEV port** for NI USB-6211 (in our example its DEV 4 but yours may differ).




Figure 19: DEV port

Next open ATIDAQFT.NET program shown below.




Figure 20: ATIDAQFT.NET

If the load cell values shown on this screen won't stop climbing then you haven't selected the proper DEV port for viewing. To change DEV ports click *Options — DAQ Hardware Options* shown below and change the device name to your DEV port. Also important to note is that our load cell works over channels 0-5 (X, Y, Z force and Rx, Ry, Rz torques).

Click *Ok* and now you should see your load cell values as shown previously. If not you have configured something wrong and you need to go back and fix this.

At this stage we have both our load cell and optical tracker setup properly and we are now ready to move on and get acquainted with our software systems for recording sensor measurements and parsing the data.




Figure 21: DAQ Hardware Options

4 Sensor Software

You have already been exposed to the native software for the load cell and optical tracker and used these native solutions to ensure our hardware is configured properly. Along with the setup software we have used thus far we will also be using three main pieces of software: return to the motion station- to record optical tracker values, loadcelly- to record load cell values, and pushpush.py to parse all the data super easily.

4.1 Return of the Motion Station

Return to motion station which was developed originally by Dr. Ferreira as part of his PhD simulator and since then has been expanded by a number of coop students. It enjoys crashing and freezing all the time but there really is no better alternative and as buggy as it is, still provides much added value.

Start by opening Return of the Motion Station. You may notice that one of the main windows doesn't enjoy life, that's fine just ignore it.




Figure 22: Return to the Motion Station

Run the project and you should be prompted with the screen shown below.

Since this is our first time running the program we should Create New however once you have already setup the experiment you can just open your already existing session folder.

Matthew Stokes




Figure 23: Create New Session Folder

Creating a new session should make the large window go black and pretty much nothing else visible will happen. (It will create all the folders structure behind the scenes among other things).




Figure 24: Return to Motion Station Running

Now we want to add our optical tracker. At this point if we havent already we can remove the stylist as we should only have the optical tracker attached to our reamer in slot 1 of our strober. I will repeat. We should only have 1 optical tracker plugged into strober in slot 1 which is the tracker attached to the reamer.

Click *Add/Configure Tracker* which will display the window shown below. We will select the *Certus* option as we are using a Certus tracker




Figure 25: Add/Configure Tracker

Selecting Certus will bring up the Certus Setup Dialog.vi shown below




Figure 26: Certus Setup

Here we want to add our smart_02 rigid body and we can name it whatever we like. This is shown in the image above. Click *Continue*.

At this stage we should be brought back to the main screen and the certus optotrak box should start beeping at us. Give this process about 2-3 minutes. Dont touch anything Motion Station is really sensitive right now. Once the process is done we should see the first slot either turn green or red which means weve successfully added our optical tracker into Motion Station.

Next we need to tell motion station to record the tracker. This is easily done with selecting the Certus 1 tracker from the list of trackers and clicking *Toggle Tracker Recording State*. This should add *ON* to the list of trackers column.

The last thing we must do is *Calibrate Tool* to add the offset we calculated earlier. Click *Calibrate Tool* to bring up the *Calibration Dialog.vi*. Insert the values you recorded into the far right column of the matrix and inset 1s along the diagonal. This is shown below.

To ensure our calibration matrix has been set, after clicking *OK* scroll down and look for the *vi Tool Calibrated Offsets* as shown in the image below. If these values arent the same as the matrix you just inputted youve done something wrong.

Wonderful! Now we are all ready to being recording. To do this as you may have guessed click the big green Start Recording button and once we stop recording we will be prompted to enter a



Figure 27: Calibration Matrix

filename. But first before we get underway with our experiment we need to configure our load cell to a state where it too is ready to record on the press of a button.

4.2 Loadcelly

Loadcelly was designed to read input voltages from a load cell on a given DEV, run these through a calibration script (given by ATI), and output values in N to screen with the ability to record these values to a file. General program can be seen below. This is the program you will use to get load cell values.




Figure 28: Loadcelly

To use this program first we must select the DEV which our load cell is on. If you dont know what DEV port this is on refer back to the step where we find the NI-6211 DEV.

Next we must load our calibration script (supplied by ATI). If you dont have this calibration script it is included as Appendix A.




Figure 29: Loadcelly Calibration File

Run the program! You should see the force and voltage values constantly jump which means yay youre reading the load cell sensor values in the program.

An option worth mentioning is the BIAS button which will bias the load cell values. You can press this if you like but my later awesome software program will handle this if you dont.

Lastly when you are ready to record you should first name the file and select a location then hit the Capture button which will turn green and start recording.

4.3 Recording

Recording is easy. Click the record button in each of the two programs then go preform your experiment. Once the experiment is done stop both recordings (motion station will prompt you to save the motion record and the file will be in XXX location)

The optical tracker file is outputted with the following headers:

Timestamp	Dx	Dy	Dz
-----------	----	----	----

Table 9: ot.csv

The load cell file is outputted with the following headers:

Matthew Stokes




Figure 30: Recording

Timestamp	Fx	Fy	Fz
-----------	----	----	----

Table 10: lc.csv

The programs can be combined but greater minds than mine are needed to do this inside of Motion Station. So instead you will have the two programs run simultaneously and you will output two text files which will then be imported into the program pushpush.py which will do many awesome things including automatically trimming the data to remove wait times before or after recording, bias, normalize, filter, and export ready to be plotted.

5 Experimental Procedure

We are looking to measure values from a push. You may record as many pushes as you like however it is recommended that you add a pause ≥ 0.25 s between each push so that the events can be considered independent.

6 Post Processing

Pushpush.py is used to extract raw information from the two text files and allow the user access to this in a much more manageable form (Motion Station outputs irrelevant information for our experiment, in a non-consistent formatting). Pushpush.py utilizes abstraction, inheritance, and the visitor design pattern to provide us easy access to our sensor data.

The main.py file is your main interaction with the program. The program should be called in the following manner and requires 3 command line arguments:

```
1 python main.py l c inputfilepath o t inputfilepath o outputfilepath
```

The output file path will be created if it does not already exist. The program will not execute if improper arguments have been provided. main.py is shown below code

You use to have to go into this main.py file and change values but that got annoying so I made a little GUI.

```
1 python main_gui.py l c inputfilepath o t inputfilepath o outputfilepath
```




Figure 31: GUI

Now it may appear as if its not working however you need to remember that most of the original timestamp values are 3,000,000,000. Therefore if we ever expect to see anything meaningful we need to normalize. Show below is an example of what the data looks like normalized.




Figure 32: GUI Normalized

Next, let's bias both the load cell and optical tracker as shown below.




Figure 33: GUI Normalized and Biased

Now this is starting to look like real data! We can remove the sections at the start and end where we have load cell values without optical tracker by clicking the *Trim* button.

Intuitively reset will reset back before we performed any operations. The last thing to touch on




Figure 34: GUI Normalized, Biased, and Trimmed

are these Excel and Excelv2.0 buttons. Excel will output the data to the output file specified in the original command. We have added some very nice functionality where we will interpolate all missing values to give each load cell and optical tracker measurements a complementary pair. This allows us to calculate stiffness and through interpolation enrich our data.

We will export to excel in the following format:

Timestamp	Fr	Dr	K
-----------	----	----	---

Table 11: Excel

The excel2.0 button will also remove all values with a Fr below 2.0N. These values can be associated with down time prepping between pushes. This allows us to quickly visually note distinctly different pushes.

All code will be included in Appendix B.

A Supplementary Software Files

A.1 smart_02

```
1 Marker Description File
2
3 Real 3D
4 3 ;Number of markers
5 1 ;Number of different views
6 Front
7 Marker      X      Y      Z      Views
8 1    -10.0955    -5.8402    0.0000    1
9 2     10.0645    -5.8402    0.0000    1
10 3     0.0309    11.6803    0.0000    1
11
12 Normals
13
14 Marker      X      Y      Z
15 0.0000    0.0000    -1.0000    1
16 0.0000    0.0000    -1.0000    1
17 0.0000    0.0000    -1.0000    1
18
19 MaxSensorError
20 0.20
21
22 Max3dError
23 0.50
24
25 MarkerAngle
26 60
27
28 3dRmsError
29 1.00
30
31 SensorRmsError
32 0.10
33
34 MinimumMarkers
35 3
36
37 MinSpread1
38 0.00
39
40 MinSpread2
41 0.00
42
43 MinSpread3
44 0.00
```

A.2 bsty0213.rig

```
1 Marker Description File
2
3 Real 3D
4 3 ;Number of markers
5 1 ;Number of different views
6 Front
7 Marker      X      Y      Z      Views
8 1      1.8639      64.4256      9.1914      1
9 2      -18.4768      63.0615      6.9789      1
10 3      -9.9560      81.3896      7.9817      1
11
12 Normals
13
14 Marker      X      Y      Z
15 -0.1079      -0.0043      0.9942      1
16 -0.1079      -0.0043      0.9942      1
17 -0.1079      -0.0043      0.9942      1
18
19 MaxSensorError
20 0.05
21
22 Max3dError
23 0.10
24
25 MarkerAngle
26 60
27
28 3dRmsError
29 0.10
30
31 SensorRmsError
32 0.05
33
34 MinimumMarkers
35 3
36
37 MinSpread1
38 0.00
39
40 MinSpread2
41 0.00
42
43 MinSpread3
44 0.00
```

A.3 FT13036Nano25ReCalibrated.cal

```
1 <?xml version="1.0" encoding="utf-8"?>
2 <!-- NOTE: To parse this file for your own software, use the "
   UserAxis" elements -->
3 <!-- to construct the calibration matrix. The "Axis" elements
   contain a scaled and -->
4 <!-- transformed version of the matrix used for ATI's internal
   purposes only. When -->
5 <!-- you send your sensor back for recalibration, or replace your
   sensor with a new -->
6 <!-- one, you will receive a new calibration file with a different
   matrix. Make -->
7 <!-- sure any custom software you write stays up to date with the
   latest version of -->
8 <!-- the calibration matrix.
   -->
9 <FTSensor Serial="FT13036" BodyStyle="Nano25" Family="DAQ" NumGages
   ="6" CalFileVersion="1.1">
10 <Calibration PartNumber="SI-125-3" CalDate="5/13/2013" ForceUnits
   ="N" TorqueUnits="N-m" DistUnits="m" OutputMode="Ground
   Referenced Differential" OutputRange="20" HWTempComp="True"
   GainMultiplier="1" CableLossDetection="False" OutputBipolar="
   True">
11 <Axis Name="Fx" values=" -0.59403 -0.44086 0.17155 32.56920
   0.24860 -34.12845 " max="125" scale="2.59826764654692"/>
12 <Axis Name="Fy" values=" 2.18077 -38.50657 -1.33492 18.50206
   -0.69315 19.67996 " max="125" scale="2.59826764654692"/>
13 <Axis Name="Fz" values=" 20.55200 -0.12973 20.36407 1.00806
   19.39320 0.18228 " max="500" scale="0.775486826825057"/>
14 <Axis Name="Tx" values=" 0.22784 -0.35271 34.87063 2.31695
   -32.82331 -0.58428 " max="3" scale="135.690306256662"/>
15 <Axis Name="Ty" values=" -40.37538 -0.02850 20.20270 0.80950
   19.10862 0.82850 " max="3" scale="135.690306256662"/>
16 <Axis Name="Tz" values=" -0.53647 -20.17550 0.95280 -19.28328
   1.35746 -19.62813 " max="3" scale="164.97900909492"/>
17 <BasicTransform Dx="0" Dy="0" Dz="0.008382" Rx="0" Ry="0" Rz="0"
   />
18 <UserAxis Name="Fx" values=" -0.22863 -0.16967 0.06602
   12.53497 0.09568 -13.13508 " max="125"/>
19 <UserAxis Name="Fy" values=" 0.83932 -14.82010 -0.51377
   7.12092 -0.26678 7.57426 " max="125"/>
20 <UserAxis Name="Fz" values=" 26.50206 -0.16729 26.25973
   1.29991 25.00777 0.23505 " max="500"/>
```

```
21     <UserAxis Name="Tx" values=" 0.00871 -0.12682 0.25268
      0.07676 -0.24413 0.05918 " max="3"/>
22     <UserAxis Name="Ty" values=" -0.29564 0.00121 0.14833
      -0.09910 0.14002 0.11620 " max="3"/>
23     <UserAxis Name="Tz" values=" -0.00325 -0.12229 0.00578
      -0.11688 0.00823 -0.11897 " max="3"/>
24     </Calibration>
25 </FTSensor>
```

B Code

B.1 Packets

B.1.1 Packet.py

```
1 class Packet(object):
2     """ Packet is an abstract class for information packets from
3         various devices"""
4
5     def __init__(self, name="Unnamed", timestamp=0):
6         """
7         Args:
8             timestamp: Timestamp from the Epoch. Default 0 (unset)
9             name: Name for this packet. Default 'Unnamed'
10
11         Attributes:
12             __name: A string containg the name of the packet
13             __timestamp: A float containing the timestamp in s from the
14                 Epoch
15         """
16         self.__name = name
17         self.__timestamp = float(timestamp)
18
19     def __str__(self):
20         return "Name: %s, Timestamp: %.5f" % (self.__name, self.
21             __timestamp)
22
23     def get_name(self): return self.__name
24     def set_name(self, name): self.__name = name
25
26     def get_timestamp(self): return self.__timestamp
27     def set_timestamp(self, value): self.__timestamp = value
```


B.1.2 OpticalTrackerData.py

```
1 from Packet import Packet
2 import math
3
4 class OpticalTrackerData(Packet):
5     """ Subclass of Packet for data captured by an optical tracker.
6     Currently holds Dx, Dy, Dz. Modify as needed
7     """
8
9     def __init__(self, name="opticaltrackerdata", timestamp=0, dx=0,
10                 dy=0, dz=0):
11         """ Super class receives timestamp value and default name of '
12             opticaltrackerdata'
13
14         Args:
15             name: packet name. opticaltrackerdata by default.
16             timestamp: Timestamp from the Epoch. Default 0 (unset)
17             dx: Displacement in X
18             dy: Displacement in Y
19             dz: Displacement in Z
20
21         Attributes:
22             __dx: A float containing displacement in X
23             __dy: A float containing displacement in Y
24             __dz: A float containing displacement in Z
25         """
26         Packet.__init__(self, name, timestamp)
27         self.__dx = float(dx)
28         self.__dy = float(dy)
29         self.__dz = float(dz)
30
31     def __str__(self):
32         return "Name: %s, Timestamp: %.5f, Dx: %.2f, Dy: %.2f, Dz: %.2f"
33             % (Packet.get_name(self), Packet.get_timestamp(self), self.
34                __dx, self.__dy, self.__dz)
35
36     def get_dx(self): return self.__dx
37     def get_dy(self): return self.__dy
38     def get_dz(self): return self.__dz
39     def get_dr(self): return math.sqrt(math.pow(self.__dx,2)+math.pow(
40         self.__dy,2)+math.pow(self.__dz,2))
41
42     def set_dx(self, value): self.__dx = value
43     def set_dy(self, value): self.__dy = value
```

```
39 | def set_dz(self, value): self.__dz = value
```

B.1.3 LoadCellData.py

```
1 from Packet import Packet
2 import math
3
4 class LoadCellData(Packet):
5     """ Subclass of Packet for data captured by a load cell.
6     Currently holds Fx, Fy, Fz. Modify as needed
7     """
8
9     def __init__(self, name="loadcelldata", timestamp=0, fx=0, fy=0, fz=0):
10         """ Super class receives timestamp value
11
12         Args:
13             name: Name of the entry. Default is "loadcelldata
14             timestamp: Timestamp from the Epoch. Default 0 (unset)
15             fx: Force in X
16             fy: Force in Y
17             fz: Force in Z
18
19         Attributes:
20             __fx: A float containing force in X
21             __fy: A float containing force in Y
22             __fz: A float containing force in Z
23         """
24         Packet.__init__(self, name, timestamp)
25         self.__fx = float(fx)
26         self.__fy = float(fy)
27         self.__fz = float(fz)
28
29     def __str__(self):
30         return "Name: %s, Timestamp: %.5f, Fx: %.2f, Fy: %.2f, Fz: %.2f"
31             % (Packet.get_name(self), Packet.get_timestamp(self), self.
32                __fx, self.__fy, self.__fz)
33
34     def get_fx(self): return self.__fx
35     def get_fy(self): return self.__fy
36     def get_fz(self): return self.__fz
37     def get_fr(self): return math.sqrt(math.pow(self.__fx,2)+math.pow(
38         self.__fy,2)+math.pow(self.__fz,2))
39
40     def set_fx(self, value): self.__fx = value
41     def set_fy(self, value): self.__fy = value
42     def set_fz(self, value): self.__fz = value
```

B.2 Sensors

B.2.1 Sensor.py

```
1 class Sensor(object):
2     """ Sensor an interface. Sensor instantiations will vary ex.
3         loadcell, optical, thermal.
4         A sensor parses data in a text file and stores data packets of
5         useful time stepped information into a list"""
6
7     def __init__(self, name="Unnamed"):
8         """
9         Args:
10             name: Name of the sensor
11
12         Attributes:
13             __name: A string containing the name of the sensor. Defaults
14                     to "Unnamed"
15             __data: A list containing packet objects for each timestamped
16                     measurement
17         """
18         self.__name = name
19         self.__data = []
20
21     def __str__(self):
22         statement = ""
23         for i in range(0, len(self.__data)):
24             statement += "Element # " + str(i) + "\n"
25             statement += str(self.__data[i]) + "\n"
26         return statement
27
28     def get_name(self): return self.__name
29     def set_name(self, name): self.__name = name
30
31     def get_data(self): return self.__data
32     def set_data(self, data_list): self.__data = data_list
33     def addData(self, packet): self.__data.append(packet)
34
35     def size(self): return len(self.__data)
```

B.2.2 LoadCell.py

```
1 from Sensor import Sensor
2 from LoadCellData import LoadCellData
3 import copy
4 class LoadCell(Sensor):
5     """A list containing all LoadCellData related to a specific
6         loadcell"""
7
8     def __init__(self, name="loadcell"):
9         Sensor.__init__(self, name)
10         self.__data = []
11
12     def __str__(self):
13         statement = ""
14         for i in range(0,len(self.__data)):
15             statement += "Element # "+str(i) + "\n"
16             statement += str(self.__data[i]) + "\n"
17         return statement
18
19     def get_data(self):
20         """ Returns a deep copy of the current data"""
21         x = copy.deepcopy(self.__data)
22         return x
23
24     def size(self): return len(self.__data)
25
26     def addData(self, packet):
27         """ Add packet to load cell. Prints an error message if packet
28             is not of
29             type LoadCellData.
30             """
31         if isinstance(packet,LoadCellData):
32             self.__data.append(packet)
33         else:
34             print "Add FAIL. NOT LOADCELLDATA packet"
35
36     def loadDataFromFile(self, filename):
37         """ Load data from an input file filename. Data is expected to
38             be in rows
39             in the format: timestamp,fx,fy,fx """
40         try:
41             reader = open(filename, 'r')
42             count = 0
43             for row in reader:
```

```
41         exploded_row = row.split(",")
42         self.__data.append(LoadCellData(filename+"_"+str(count),
43             exploded_row[0], exploded_row[1], exploded_row[2],
44             exploded_row[3]))
45         count = count + 1
46     except:
47         print "FAIL loading load cell data from file"
48
49 def accept(self, data_exchange):
50     """ Visitor design pattern. Allow for a visit. """
51     data_exchange.visit(loadcell=self)
```

B.2.3 OpticalTracker.py

```
1 from Sensor import Sensor
2 from OpticalTrackerData import OpticalTrackerData
3 import copy
4
5 class OpticalTracker(Sensor):
6     """ Instantiation of Sensor. Stores a list of OpticalTrackerData.
7     """
8
9     def __init__(self, name="Optical Tracker"):
10         """ Initializes "Optical Tracker" as name for every new
11         OpticalTracker Sensor
12
13         Attributes:
14             __data: List of OpticalTrackerData
15         """
16         Sensor.__init__(self, name)
17         self.__data = []
18
19     def __str__(self):
20         statement = ""
21         for i in range(0, len(self.__data)):
22             statement += "Element # " + str(i) + "\n"
23             statement += str(self.__data[i]) + "\n"
24         return statement
25
26     def get_data(self):
27         """get a deep copy of the current data list"""
28         return copy.deepcopy(self.__data)
29
30     def size(self): return len(self.__data)
31
32     def addData(self, packet):
33         """ Add a new OpticalTrackerData to the list
34         Args:
35             packet: OpticalTrackerData to add to the list
36         Return:
37             Prints message in error
38         """
39         if isinstance(packet, OpticalTrackerData):
40             self.__data.append(packet)
41         else:
42             print "Add FAIL. Packet NOT OF TYPE OpticalTrackerData"
```

```

42 def loadDataFromFile(self, filename):
43     """ Loads data from filename into OpticalTrackerData objects
        than stores these objects in the OpticalTracker list.
        Extremely odd format outputted by LABVIEW
44     Args:
45         loadcell_file: Input file for loadcell
46         optical_tracker_file: Input file for optical tracker
47     """
48     try:
49         reader = open(filename, 'r')
50         reader.readline() # pass first line
51         count = 0
52         for row in reader:
53             exploded_row = row.split(" ")
54             self.__data.append(OpticalTrackerData(filename+"_"+str(count)
                ),exploded_row[0].split('\t')[0],exploded_row[3],
                exploded_row[7],exploded_row[11]))
55             count = count + 1
56     except:
57         print "FAIL loading optrical tracker data from file"
58
59 def accept(self, data_exchange):
60     """ Visitor design pattern. Allow visit. """
61     data_exchange.visit(opticaltracker=self)

```


B.3 Visitor Design Pattern

B.3.1 Reamer.py

```
1 from OpticalTracker import OpticalTracker
2 from LoadCell import LoadCell
3 from DataExchange import DataExchange
4
5 class Reamer(object):
6     """ A Reamer is a collection of sensors. """
7
8     def __init__(self, opticaltracker=OpticalTracker(), loadcell=
9         LoadCell()):
10         """
11         Args:
12             opticaltracker: An optical tracker. Will create empty optical
13                             tracker by default
14             loadcell: A load cell. Will create empty load cell by default
15
16         Attributes:
17             __ot: A list of optical trackers
18             __lc: A list of load cells
19         """
20         if not isinstance(opticaltracker, OpticalTracker) or not
21             isinstance(loadcell, LoadCell):
22             print "Improper initialization of Reamer. Should be types
23                 OpticalTracker and LoadCell"
24         else:
25             self.__ot = opticaltracker
26             self.__lc = loadcell
27
28     def __str__(self):
29         return str(self.__ot) + str(self.__lc)
30
31     def loadDataFromFile(self, loadcell_file, optical_tracker_file):
32         """ Loads data from input file into loadcell and opticaltracker
33             lists respectfully. Requires one load cell and one optical
34             tracker file
35
36         Args:
37             loadcell_file: Input file for loadcell
38             optical_tracker_file: Input file for optical tracker
39         """
40         self.__ot.loadDataFromFile(optical_tracker_file)
41         self.__lc.loadDataFromFile(loadcell_file)
42
43     def accept(self, data_exchange):
```

```
37     """ Visitor design pattern. Allows for a DataExchange object to
38         visit to extract information
39     Args:
39         data_exchange: DataExchange object. Called via: Reamer.visit(
40             DataExchange Obj)
40     """
41     data_exchange.visit(loadcell=self.__lc, opticaltracker=self.__ot
42         )
```

B.3.2 DataExchange.py

```
1 import sys
2 from LoadCell import LoadCell
3
4 class DataExchange(object):
5     """ Abstact visitor class """
6     def __init__(self):
7         pass
8
9     def visit(self, sensor):
10         pass
```

B.3.3 Chart.py

```
1 import sys
2 import math
3 import os
4 import copy
5 from LoadCell import LoadCell
6 from DataExchange import DataExchange
7 from LoadCellData import LoadCellData
8 from OpticalTrackerData import OpticalTrackerData
9 from OpticalTracker import OpticalTracker
10 from StiffnessElement import StiffnessElement
11
12 #from pylab import *
13 #import matplotlib.pyplot as plt
14
15 class Chart(DataExchange):
16     def __init__(self):
17         self.__lc = [] # list of lists
18         self.__ot = [] # list of lists
19         self.__x = []
20         self.__k = []
21
22     def visit(self, loadcell=None, opticaltracker=None):
23         if loadcell != None:
24             self.addLoadCell(loadcell.get_data())
25         if opticaltracker != None:
26             self.addOpticalTracker(opticaltracker.get_data())
27
28     def addLoadCell(self, loadcell):
29         self.__lc.append(copy.deepcopy(loadcell))
30         for i in range(0, len(loadcell)):
31             self.__x.append(copy.deepcopy(loadcell[i]))
32         self.__x = sorted(self.__x, key=lambda sensor: sensor.get_timestamp())
33
34     def addOpticalTracker(self, ot):
35         self.__ot.append(copy.deepcopy(ot))
36         for i in range(0, len(ot)):
37             self.__x.append(copy.deepcopy(ot[i]))
38         self.__x = sorted(self.__x, key=lambda sensor: sensor.get_timestamp())
39
40     def normalize(self):
41         start = float(self.__x[0].get_timestamp())
```

```

42     for i in range(0,len(self.__x)):
43         self.__x[i].set_timestamp(float(self.__x[i].get_timestamp()) -
44                                     start)
45
46 def printTimestamp(self):
47     for item in self.__x:
48         print item.get_timestamp()
49
50 def reset(self):
51     self.__x = []
52
53     # add from _ot and _lc
54     for loadcell in self.__lc:
55         for i in range(0,len(loadcell)):
56             self.__x.append(copy.deepcopy(loadcell[i]))
57     for opticaltracker in self.__ot:
58         for i in range(0,len(opticaltracker)):
59             self.__x.append(copy.deepcopy(opticaltracker[i]))
60
61     self.__x = sorted(self.__x, key=lambda sensor: sensor.
62                       get_timestamp())
63
64 def printLC(self):
65     for item in self.__x:
66         if isinstance(item, LoadCellData):
67             print item.get_fx()
68             print item.get_fy()
69             print item.get_fz()
70
71 def printOT(self):
72     for item in self.__x:
73         if isinstance(item, OpticalTrackerData):
74             print item.get_dr()
75
76 def selectRange(self, start, end):
77     """ start and end in seconds """
78     if start > end:
79         print "ERROR: Select Range - Starting value is greater than
80             ending value"
81         return
82     start_value = float(self.__x[0].get_timestamp())
83     end_value = float(self.__x[-1].get_timestamp())

```

```

84     start_counter = 0
85     end_counter = 0
86
87     for packet in self.__x:
88         if packet.get_timestamp() <= start + start_value:
89             start_counter = start_counter + 1
90         if packet.get_timestamp() <= end + start_value:
91             end_counter = end_counter + 1
92
93     self.__x = self.__x[start_counter:end_counter]
94
95     def duration(self):
96         return self.__x[-1].get_timestamp() - self.__x[0].get_timestamp()
97
98     def trimStart(self, seconds):
99         start = float(self.__x[0].get_timestamp())
100
101         position = 0
102         for packets in self.__x:
103             if start + seconds < packets.get_timestamp():
104                 break
105             else:
106                 position = position + 1
107
108         self.__x = self.__x[position:]
109
110     def trimEnd(self, seconds):
111         end = float(self.__x[-1].get_timestamp())
112
113         position = 0
114         for packets in self.__x:
115             if end - seconds < packets.get_timestamp():
116                 break
117             else:
118                 position = position + 1
119
120         self.__x = self.__x[0:position]
121
122     def autoTrim(self):
123
124         # start point
125         start = 1
126         if isinstance(self.__x[0], LoadCellData):
127             while isinstance(self.__x[start], LoadCellData):

```

```

128         start = start + 1
129     elif isinstance(self.__x[0], OpticalTrackerData):
130         while isinstance(self.__x[start], OpticalTrackerData):
131             start = start + 1
132
133     end = len(self.__x) - 1
134     if isinstance(self.__x[end], LoadCellData):
135         while isinstance(self.__x[end], LoadCellData):
136             end = end - 1
137     elif isinstance(self.__x[end], OpticalTrackerData):
138         while isinstance(self.__x[end], OpticalTrackerData):
139             end = end - 1
140
141     self.__x = self.__x[start-1:end+2]
142
143     def biasLoadCell(self):
144         start_counter = 0
145         while not isinstance(self.__x[start_counter], LoadCellData):
146             start_counter = start_counter + 1
147
148         start_x = float(self.__x[start_counter].get_fx())
149         start_y = float(self.__x[start_counter].get_fy())
150         start_z = float(self.__x[start_counter].get_fz())
151
152         for packets in self.__x:
153             if isinstance(packets, LoadCellData):
154                 packets.set_fx(packets.get_fx() - start_x)
155                 packets.set_fy(packets.get_fy() - start_y)
156                 packets.set_fz(packets.get_fz() - start_z)
157
158     def biasOpticalTracker(self):
159         start_counter = 0
160         while not isinstance(self.__x[start_counter], OpticalTrackerData):
161             start_counter = start_counter + 1
162
163         start_x = float(self.__x[start_counter].get_dx())
164         start_y = float(self.__x[start_counter].get_dy())
165         start_z = float(self.__x[start_counter].get_dz())
166
167         for packets in self.__x:
168             if isinstance(packets, OpticalTrackerData):
169                 packets.set_dx(packets.get_dx() - start_x)
170                 packets.set_dy(packets.get_dy() - start_y)
171                 packets.set_dz(packets.get_dz() - start_z)

```

```

172
173 def printX(self):
174     for item in self.__x:
175         print item
176
177 def printOpticalTrackingData(self):
178     for sensor in self.__ot:
179         for packets in sensor:
180             print packets
181
182 def toExcel(self, filename):
183     try:
184
185         if not os.path.exists(os.path.abspath(os.path.join(filename,
186             os.pardir)))):
187             os.makedirs(os.path.abspath(os.path.join(filename, os.pardir
188                 )))
189
190         f = open(filename, 'w')
191         f.write("Timestamp,Resultant_Force,Resultant_Displacement\n")
192         for packet in self.__x:
193             line = str(packet.get_timestamp()) + ","
194             if isinstance(packet, LoadCellData):
195                 line += str(packet.get_fr())
196             elif isinstance(packet, OpticalTrackerData):
197                 line += "," + str(packet.get_dr())
198             f.write(line+"\n")
199         except IOError as e:
200             print "Error opening output file for writing...impressive"
201
202 def calculateStiffness(self):
203     self.__interpolate()
204
205 def isolate1push(self, threshold):
206     if self.__k == []:
207         self.__interpolate()
208
209     new_k = []
210
211     for element in self.__k:
212         if element.get_fr() > threshold:
213             new_k.append(element)
214
215     self.__k = new_k

```



```

215 def toExcelInterpolated(self, filename):
216     try:
217
218         if not os.path.exists(os.path.abspath(os.path.join(filename,
219             os.pardir)))):
220             os.makedirs(os.path.abspath(os.path.join(filename, os.pardir
221                 )))
222
223         f = open(filename, 'w')
224         f.write("Timestamp,Resultant_Force,Resultant_Displacement,
225             Stiffness\n")
226         if self.__k == []:
227             self.__interpolate()
228
229         for element in self.__k:
230             line = str(element.get_timestamp()) + ","
231             line += str(element.get_fr()) + ","
232             line += str(element.get_dr()) + ","
233             line += str(element.get_stiffness())
234             f.write(line+"\n")
235
236         print "toExcelInterpolated successful"
237
238     except IOError as e:
239         print "Error opening output file for writing...impressive"
240
241 def toExcelPushSeperated(self, filename):
242     try:
243         if not os.path.exists(os.path.abspath(os.path.join(filename,
244             os.pardir)))):
245             os.makedirs(os.path.abspath(os.path.join(filename, os.pardir
246                 )))
247
248         f = open(filename, 'w')
249         f.write("Timestamp,Resultant_Force,Resultant_Displacement,
250             Stiffness\n")
251         if self.__k == []:
252             self.__interpolate()
253
254         prev = self.__k[0].get_timestamp()
255         for element in self.__k:
256             line = str(element.get_timestamp()) + ","
257             line += str(element.get_fr()) + ","
258             line += str(element.get_dr()) + ","
259             line += str(element.get_stiffness())

```

```

254         if element.get_timestamp() - prev > .25:
255             for i in range(0,5): f.write("\n")
256             f.write(line+"\n")
257             prev = element.get_timestamp()
258
259         print "toExcelPushSeperated successful"
260
261     except IOError as e:
262         print "Error opening output file for writing...impressive"
263
264
265     def dumpTimeLC(self):
266         x = []
267         for item in self.__x:
268             if isinstance(item, LoadCellData):
269                 x.append(item.get_timestamp())
270         return x
271
272     def dumpTimeOT(self):
273         x = []
274         for item in self.__x:
275             if isinstance(item, OpticalTrackerData):
276                 x.append(item.get_timestamp())
277         return x
278
279     def dumpLC(self):
280         x = []
281         for item in self.__x:
282             if isinstance(item, LoadCellData):
283                 x.append(item.get_fr())
284         return x
285
286     def dumpOT(self):
287         x = []
288         for item in self.__x:
289             if isinstance(item, OpticalTrackerData):
290                 x.append(item.get_dr())
291         return x
292
293     def isolate(self, start, end):
294         self.__x = self.__x[start:end]
295
296     ''' 2 pass: 1) loadcell interp 2) optical tracker interp
297         returns a list of turpals'''
298     def __interpolate(self):

```

```

299
300     stiff_list = []
301
302     # pass 1: loadcell interpolation
303     interp1 = False # bool
304     interp1_val = 0
305     interp2 = False
306     interp2_val = 0
307
308     for loc in range(0, len(self.__x)):
309
310         # if found a load cell value
311         if isinstance(self.__x[loc], LoadCellData):
312             if interp2 == True:
313                 interp2_val = loc
314                 self.__interpolateLC(interp1_val, interp2_val, stiff_list)
315                 interp2 = False
316
317             if interp2 == False:
318                 interp1 = True
319                 interp1_val = loc
320                 stiff_list.append(StiffnessElement(timestamp=self.__x[loc]
321                     ].get_timestamp(), fr=self.__x[loc].get_fr())) # add lc
322                     data to stiffness element
321         else: # optical tracker element
322             interp2 = True
323
324
325     # pass 2: optical tracker interpolation
326     interp1 = False
327     interp1_val = 0
328     interp2 = False
329     interp2_val = 0
330
331     for loc in range(0, len(self.__x)):
332
333         # if found an optical tracker value
334         if isinstance(self.__x[loc], OpticalTrackerData):
335             if interp2 == True:
336                 interp2_val = loc
337                 self.__interpolateOT(interp1_val, interp2_val, stiff_list)
338                 interp2 = False
339
340             if interp2 == False:
341                 interp1 = True

```

```

342         interp1_val = loc
343     try:
344         stiff_list[loc].set_dr(self.__x[loc].get_dr()) # add lc
345             data to stiffness element
346     except IndexError:
347         pass
348
349     else: # optical tracker element
350         interp2 = True
351
352     self.__k = stiff_list
353
354     ''' interpolate values using timestamp '''
355     def __interpolateLC(self, start, end, stiff_list):
356         for element in range(start+1, end):
357
358             # y = y_a + (y_b - y_a)[(x-x_a)/(x_b-x_a)]
359             interp_fr = self.__x[start].get_fr() + ((self.__x[end].get_fr
360                 () - self.__x[start].get_fr()) * (( self.__x[element].
361                 get_timestamp() - self.__x[start].get_timestamp()) / (self.
362                 __x[end].get_timestamp() - self.__x[start].get_timestamp())
363             ))
364             #print "Element: " + str(element) + " Value: " + str(interp_fr
365             )
366
367             stiff_list.append(StiffnessElement(timestamp=self.__x[element
368                 ].get_timestamp(), fr=interp_fr))
369
370     ''' interpolate values using timestamp '''
371     def __interpolateOT(self, start, end, stiff_list):
372         for element in range(start+1, end):
373
374             # y = y_a + (y_b - y_a)[(x-x_a)/(x_b-x_a)]
375             interp_dr = self.__x[start].get_dr() + ((self.__x[end].get_dr
376                 () - self.__x[start].get_dr()) * (( self.__x[element].
377                 get_timestamp() - self.__x[start].get_timestamp()) / (self.
378                 __x[end].get_timestamp() - self.__x[start].get_timestamp())
379             ))
380             #print "Element OT: " + str(element) + " Value: " + str(
381                 interp_dr)
382             try:
383                 stiff_list[element].set_dr(interp_dr)
384             except IndexError:
385                 pass # last element value is undetermined

```

B.3.4 StiffnessElement.py

```
1
2 class StiffnessElement(object):
3
4     def __init__(self, timestamp=0, fr=0, dr=0):
5         self.__timestamp = timestamp
6         self.__fr = fr
7         self.__dr = dr
8
9     def __str__(self):
10         return "Timestamp: %s, Fr: %.2f, Dr: %.2f" % (self.__timestamp,
11             self.__fr, self.__dr)
12
13     def set_fr(self, value):
14         self.__fr = value
15
16     def set_dr(self, value):
17         self.__dr = value
18
19     def set_timestamp(self, value):
20         self.__timestamp = value
21
22     def get_timestamp(self): return self.__timestamp
23     def get_fr(self): return self.__fr
24     def get_dr(self): return self.__dr
25
26     def get_stiffness(self):
27         if self.__dr!=0:
28             return (self.__fr/self.__dr) * 1000 # N/m
29         else:
30             return 0
```

B.4 GUI

B.4.1 gui.py

```
1 import numpy as np
2 import matplotlib.pyplot as plt
3 from matplotlib.widgets import Button
4
5 class Index:
6
7     # takes a chart object and preforms necessary actions
8     def __init__(self, chart, ax, output_file):
9         self.__chart = chart
10        self.__ax = ax
11        self.__output_file = output_file
12        self.__plot_init()
13
14        # data for line values
15        x = self.__chart.dumpTimeLC()
16        y = self.__chart.dumpLC()
17
18        x2 = self.__chart.dumpTimeOT()
19        y2 = self.__chart.dumpOT()
20
21        # plotting data
22        self.__l1, = self.__ax.plot(x, y, lw=2)
23        self.__l2, = self.__ax.plot(x2, y2, lw=2)
24
25        self.__isolated = 0
26
27    def __reload(self):
28
29        self.__l1.set_xdata(self.__chart.dumpTimeLC())
30        self.__l1.set_ydata(self.__chart.dumpLC())
31
32        self.__l2.set_xdata(self.__chart.dumpTimeOT())
33        self.__l2.set_ydata(self.__chart.dumpOT())
34
35        self.__ax.relim()
36        self.__ax.autoscale()
37
38        plt.draw()
39
40    def autoTrim(self, event):
41        print "Autotrim"
42        self.__chart.autoTrim()
```

```

43     self.__reload()
44
45 def normalize(self, event):
46     print "Normalize"
47     self.__chart.normalize()
48     self.__reload()
49
50 def biasLoadCell(self, event):
51     print "Bias Load Cell"
52     self.__chart.biasLoadCell()
53     self.__reload()
54
55 def biasOpticalTracker(self, event):
56     print "Bias Optical Tracker"
57     self.__chart.biasOpticalTracker()
58     self.__reload()
59
60 def isolate(self, event):
61     print "Isolate"
62     self.__chart.isolate1push(4)
63     self.__reload()
64
65 def toExcelInterpolated(self, event):
66     print "to Excel Interpolated"
67     self.__chart.toExcelInterpolated(self.__output_file)
68     self.__reload()
69
70 def toExcelPushSeperated(self, event):
71     print "to Excel Push Seperated"
72     self.__chart.toExcelPushSeperated(self.__output_file)
73     self.__reload()
74
75 def reset(self, event):
76     print "Reset"
77     self.__chart.reset()
78     self.__reload()
79
80 def __plot_init(self, fontsize=12):
81     self.__ax.plot([1])
82     self.__ax.locator_params(nbins=3)
83     self.__ax.set_xlabel('Time', fontsize=fontsize)
84     self.__ax.set_ylabel('Force or Displacement', fontsize=fontsize)
85     self.__ax.set_title('Reamer Results', fontsize=fontsize)
86
87 def isolate(self, val):

```

```
88     print int(val)
89     if int(val) != self.__isolated:
90         self.__chart.isolate1push(val)
91         print "change"
92         self.__isolated = int(val)
93
94     self.__reload()
```


B.4.2 main_gui.py

```
1 from Reamer import Reamer
2 from Chart import Chart
3
4 import argparse
5 import sys
6
7 from gui import Index
8 import numpy as np
9 import matplotlib.pyplot as plt
10 from matplotlib.widgets import Button, Slider
11
12
13 ##### Command Line Arguments #####
14
15 # get command line arguments input, output, loadcel
16 parser = argparse.ArgumentParser()
17 parser.add_argument("-lc", "--loadcell", help="file containing load
    cell data")
18 parser.add_argument("-ot", "--opticaltracker", help="file containing
    optical tracker data")
19 parser.add_argument("-o", "--output", help="output file")
20 args = parser.parse_args()
21
22 # error check command line arguments
23 if not args.loadcell:
24     parser.error("please include loadcell input file")
25 if not args.opticaltracker:
26     parser.error("please include optical tracker input file")
27 if not args.output:
28     parser.error("pleasae include output file")
29
30
31 ##### Visitor Pattern Setup #####
32
33 # Create objects reamer and chart
34 r = Reamer()
35 c = Chart()
36
37
38 # load data from file(s) into reamer
39 r.loadDataFromFile(args.loadcell, args.opticaltracker)
40
41 # setup visitor
```

```

42 r.accept(c)
43
44
45 ##### GUI Setup #####
46
47
48 fig = plt.figure()
49
50 ax1 = plt.subplot2grid((5, 5), (0, 0), colspan=5, rowspan=3)
51 pos_trim = plt.subplot2grid((5, 5), (3, 0), colspan=1, rowspan=1)
52 pos_normalize = plt.subplot2grid((5, 5), (3, 1), colspan=1, rowspan
    =1)
53 pos_bias_lc = plt.subplot2grid((5, 5), (3, 2), colspan=1, rowspan=1)
54 pos_bias_ot = plt.subplot2grid((5, 5), (3, 3), colspan=1, rowspan=1)
55 #pos_isolate = plt.subplot2grid((5, 5), (3, 4), colspan=1, rowspan
    =1)
56 pos_excel_interp = plt.subplot2grid((5, 5), (4, 0), colspan=1,
    rowspan=1)
57 pos_excel_seperated = plt.subplot2grid((5, 5), (4, 1), colspan=1,
    rowspan=1)
58 pos_reset = plt.subplot2grid((5, 5), (3, 4), colspan=1, rowspan=1)
59
60 plt.tight_layout()
61
62 # button names
63 b_trim = Button(pos_trim, 'Trim', hovercolor='0.25')
64 b_normalize = Button(pos_normalize, 'Normalize', hovercolor='0.25')
65 b_bias_lc = Button(pos_bias_lc, 'Bias LC', hovercolor='0.25')
66 b_bias_ot = Button(pos_bias_ot, 'Bias OT', hovercolor='0.25')
67 b_excel_interp = Button(pos_excel_interp, 'Excel', hovercolor='0.25'
    )
68 b_excel_seperated = Button(pos_excel_seperated, 'Excelv2.0',
    hovercolor='0.25')
69 b_reset = Button(pos_reset, 'Reset', hovercolor='0.25')
70
71 #sl_isolate = Slider(pos_isolate, 'Isolate', valmin=0, valmax=5,
    valinit=0, valfmt='%d')
72
73
74 # button on_click
75 callback = Index(c, ax1, args.output)
76
77 b_trim.on_clicked(callback.autoTrim)
78 b_normalize.on_clicked(callback.normalize)
79 b_bias_lc.on_clicked(callback.biasLoadCell)

```

```
80 b_bias_ot.on_clicked(callback.biasOpticalTracker)
81 b_excel_interp.on_clicked(callback.toExcelInterpolated)
82 b_excel_seperated.on_clicked(callback.toExcelPushSeperated)
83 b_reset.on_clicked(callback.reset)
84 #sl_isolate.on_changed(callback.isolate)
85
86 ##### Running GUI #####
87
88 plt.show()
```